




GTS

GUIDED THERAPY SYSTEM

New intuitive User Interface

Unique GTS: fine tune the treatment to your patient's condition

The intelligent Shockmaster patient Guided Therapy System will automatically adapt treatment level based on your patient's condition.



New Sub Menu's
based on affected tissue

Application Screen

New Design
Look & Feel in line with
the new 400 series

Treatment Screen
Plug & Play

Contra-indications

included in the
ShockMaster 300.

Pop up for revision

Appears after 1.25 million shocks
Timely revision is essential for a reliable and
effective functioning of the ShockMaster device.

Safety restrictions

to prevent wrong
handpiece usage.



New pathologies

- Spasticity in children after Cerebral Palsy
- Trigger Finger
- Scar Tissue
- Biceps Tendinopathy
- Erectile Dysfunction



Backwards compatibility

1 Software

Easy software update via USB port is
available for all ShockMaster devices.

2 Handpiece

The short handpiece may be used on all
ShockMaster devices provided that the
mandatory software update is performed.